

Product Datasheet

Testosterone-3-HRP Antigen 170068

Product Name	Testosterone-3-HRP Antigen
Catalog Number	170068
Description	Synthetic Testosterone-3, HRP conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Storage	Storage at 2-8 °C. Do not freeze.
Note	This product was previously available with catalog code P91-92-71H. It has been produced using the same conjugation protocol.
Buffer	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016% Gentamicin sulfate, 0.005% thimerosal
Shelf Life	Unspecified