

## Product Datasheet

### Testosterone-3-HRP Antigen 170068

<b>Product Name</b>	Testosterone-3-HRP Antigen
<b>Catalog Number</b>	170068
<b>Description</b>	Synthetic Testosterone-3, HRP conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Storage at 2-8 °C. Do not freeze.
<b>Note</b>	This product was previously available with catalog code P91-92-71H. It has been produced using the same conjugation protocol.
<b>Buffer</b>	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016% Gentamicin sulfate, 0.005% thimerosal
<b>Shelf Life</b>	Unspecified