

Product Datasheet

Anti-Human CD130 Monoclonal Antibody, Azide Free Clone B-T2 853.860

Product Name	Anti-Human CD130 Monoclonal Antibody, Azide Free Clone B-T2
Catalog Number	853.860
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	Natural soluble CD130 (gp130)
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises Gp130, common subunit for IL-6, IL-11, OSM, LIF, CNTF, CT-1 receptors, a 130-140 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells