

## Product Datasheet

### Amphetamine-HRP Antigen 170100

<b>Product Name</b>	Amphetamine-HRP Antigen
<b>Catalog Number</b>	170100
<b>Description</b>	Synthetic Amphetamine, HRP conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid; Antigen-HRP conjugate
<b>Storage</b>	2-8 °C. Do not freeze
<b>Note</b>	This product was previously available with Catalog code P91-99-71H. It has been produced using the same conjugation protocol. Linked via para position
<b>Associated Products</b>	C03-99-72M-P
<b>Buffer</b>	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016 % Gentamicin sulfate, 0.005 % thimerosal
<b>Matched Pairs</b>	C03-99-72M-P (capture)
<b>Shelf Life</b>	Unspecified, storage at 2-8 °C. Do not freeze