

Product Datasheet

Phencyclidine (PCP)-HRP Antigen 170090-1

Product Name	Phencyclidine (PCP)-HRP Antigen
Catalog Number	170090
Description	Synthetic Phencyclidine, HRP conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-HRP conjugate
Storage	2-8 °C. Do not freeze
Note	This product was previously available with Catalog code P91-99-61H. It has been produced using the same conjugation protocol.
Associated Products	C03-99-72M-P
Buffer	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016 % Gentamicin sulfate, 0.005 % thimerosal
Shelf Life	Unspecified, storage at 2-8 °C. Do not freeze