

Product Datasheet

Serum, Direct Bilirubin Samples 0.1 - 0.4 mg/dL mg/dL 991-24-PS-DBIL2

Product Name Serum, Direct Bilirubin Samples 0.1 - 0.4 mg/dL mg/dL

Catalog Number 991-24-PS-DBIL2

Brand Lee Biosolutions

Form/Appearance Liquid

Storage -20°C

Source Single Human Donor

Shelf Life Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.