

## Product Datasheet

### High Density Lipoprotein (HDL) Serum 991-24-PS-HDL

<b>Product Name</b>	High Density Lipoprotein (HDL) Serum
<b>Catalog Number</b>	991-24-PS-HDL
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20°C
<b>Source</b>	Single Human Donor
<b>Shelf Life</b>	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.