

Product Datasheet

Rheumatoid Factor (RF) from Human Plasma; ? 25,000 IU/mL 508-15

Product Name	Rheumatoid Factor (RF) from Human Plasma; \geq 25,000 IU/mL
Catalog Number	508-15
Tested Applications	Control, Calibration, QC Materials
Brand	Lee Biosolutions
Form/Appearance	Gray to brown, hazy liquid
Concentration	\geq 25,000 IU/mL (ELISA)
Storage	-20°C
Note	Protein: Reported in mg/mL (Biuret) Specific Activity: Reported in IU/mg protein
Source	Human Plasma
Buffer	20 mM sodium phosphate, 0.5 M sodium chloride, pH 7.2
Shelf Life	1 year from date of manufacture
Preservatives	0.02 % ProClin 300