

## Product Datasheet

### Rheumatoid Factor (RF) from Human Plasma; ? 25,000 IU/mL 508-15

<b>Product Name</b>	Rheumatoid Factor (RF) from Human Plasma; $\geq$ 25,000 IU/mL
<b>Catalog Number</b>	508-15
<b>Tested Applications</b>	Control, Calibration, QC Materials
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Gray to brown, hazy liquid
<b>Concentration</b>	$\geq$ 25,000 IU/mL (ELISA)
<b>Storage</b>	-20°C
<b>Note</b>	Protein: Reported in mg/mL (Biuret) Specific Activity: Reported in IU/mg protein
<b>Source</b>	Human Plasma
<b>Buffer</b>	20 mM sodium phosphate, 0.5 M sodium chloride, pH 7.2
<b>Shelf Life</b>	1 year from date of manufacture
<b>Preservatives</b>	0.02 % ProClin 300