

Product Datasheet

Anti-Cryptosporidium V7641

Product Name	Anti-Cryptosporidium
Catalog Number	V7641
Species/Host	Rabbit
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from culture medium by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Note	Periodate sensitive
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Giardia lamblia, Candida albicans & human cells
Specificity	Cyst cell wall glyco-antigen
Matched Pairs	ELISA pairs (capture-conjugate) : V7641-V7642, V7642-V7641
Shelf Life	Four years from date of manufacture