

## Product Datasheet

### Anti-Cryptosporidium V7642

<b>Product Name</b>	Anti-Cryptosporidium
<b>Catalog Number</b>	V7642
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA, IFA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from culture medium by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Note</b>	Periodate sensitive
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Non-reactive with Giardia lamblia, Candida albicans & human cells
<b>Specificity</b>	Cyst cell wall glyco-antigen
<b>Matched Pairs</b>	ELISA pairs (capture-conjugate) : V7641-V7642, V7642-V7642, V7642-V7641
<b>Shelf Life</b>	Four years from date of manufacture