

Product Datasheet

Anti-Fungal beta-glucan V8201

Product Name	Anti-Fungal beta-glucan
Catalog Number	V8201
Species/Host	Rabbit
Brand	ViroStat
Form/Appearance	Purified from culture medium by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Note	Extraction of specimen with 0.4M NaOH required for optimum sensitivity.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide
Cross Reactivity	Reactive with Candida (4 species), Aspergillus (4 species), Penicillium notatum, Saccharomyces cerevisiae, Pichia pastoris, and Mucor hiemalis. Others not tested.
Specificity	Numerous fungal species
Shelf Life	Four years from date of manufacture