

Product Datasheet

Anti-Streptococcus pneumoniae V0494

Product Name Anti-Streptococcus pneumoniae

Catalog Number V0494

Species/Host Rabbit

Tested Applications LF, ELISA, IFA

Brand ViroStat

Form/Appearance Purified from culture medium by protein A chromatography

Concentration ≥ 1 mg/ml

Storage 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.

Isotype IgG1

Clonality Monoclonal

Purity $\geq 90\%$

Buffer 0.01M Phosphate Buffered Saline, pH 7.5, containing 0.1% sodium azide

Cross Reactivity Not determined

Specificity Common cell wall polysaccharide (urinary antigen) of *S. pneumoniae*

Matched Pairs ELISA pairs (capture-conjugate): V0494-V0495

Shelf Life Four years from date of manufacture