

Product Datasheet

Anti-Canine Heartworm V3221

Product Name	Anti-Canine Heartworm
Catalog Number	V3221
Species/Host	Goat
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is not absorbed for related filarial Nematodes. Non-reactive with canine serum
Specificity	Soluble & structural antigens of the adult worm
Shelf Life	Six years from date of manufacture