

Product Datasheet

Anti-Papillomavirus V5001

Product Name	Anti-Papillomavirus
Catalog Number	V5001
Species/Host	Goat
Tested Applications	ELISA, IHC, IFA
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HPV types 16 & 11 (other types not tested). Reactive in ELISA with Bovine papillomavirus & HPV VLPS. Non-reactive with normal tissue (IHC) & cultured cells (IFA)
Specificity	group
Shelf Life	Six years from date of manufacture