

Product Datasheet

Anti-Herpes simplex virus, types 1&2 HRP V0104

Product Name	Anti-Herpes simplex virus, types 1&2 HRP
Catalog Number	V0104
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Reactive with types 1 & 2. Non-reactive with HEp-2 cells
Specificity	ICPs + late structural (virion) antigens
Shelf Life	Two years from date of manufacture