

Product Datasheet

Anti-Epstein-Barr virus (VCA-p18) V0291

Product Name	Anti-Epstein-Barr virus (VCA-p18)
Catalog Number	V0291
Species/Host	Goat
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	p18 of VCA
Shelf Life	Six years from date of manufacture