

Product Datasheet

Anti-Leptospira V0401

Product Name	Anti-Leptospira
Catalog Number	V0401
Species/Host	Rabbit
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive in ELISA with <i>L. biflexa</i> , <i>L. australis</i> , <i>L. canicola</i> , <i>L. grippityphosa</i> , <i>L. icterohaemorrhagiae</i> , <i>L. pomona</i> , <i>L. serjroe</i> . Antiserum has not been absorbed and may react with related microorganisms
Specificity	All antigens
Shelf Life	Six years from date of manufacture