

Product Datasheet

Anti-E. coli V1001

| | |
|----------------------------|--|
| Product Name | Anti-E. coli |
| Catalog Number | V1001 |
| Species/Host | Rabbit |
| Tested Applications | ELISA |
| Brand | ViroStat |
| Form/Appearance | Purified from antiserum by protein A chromatography |
| Concentration | 5 mg/ml |
| Storage | Short term (<6 months): 2-8°C Long term: -20°C |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Antiserum is not absorbed and does react with related Enterobacteriaceae |
| Specificity | All "O" and "K" antigenic serotypes of Escherichia coli |
| Shelf Life | Six years from date of manufacture |