

Product Datasheet

Anti-Influenza B virus V1101

| | |
|----------------------------|--|
| Product Name | Anti-Influenza B virus |
| Catalog Number | V1101 |
| Species/Host | Goat |
| Tested Applications | ELISA, WB, IFA, Indirect immunofluorescence |
| Brand | ViroStat |
| Form/Appearance | Purified IgG fraction prepared by methods developed by Medix Biochemica |
| Concentration | 5 mg/ml |
| Storage | Short term (<6 months): 2-8°C Long term: -20°C |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Non-reactive with Influenza A, Para 1-3, RSV & Adenovirus. Non-reactive with monkey kidney cells by indirect immunofluorescence |
| Specificity | Purified virions |
| Shelf Life | Six years from date of manufacture |