

Product Datasheet

Anti-Influenza B virus Biotin V1107

Product Name	Anti-Influenza B virus Biotin
Catalog Number	V1107
Species/Host	Goat
Tested Applications	ELISA, WB, IFA, Indirect immunofluorescence
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Influenza A, Para 1-3, RSV & Adenovirus. Non-reactive with monkey kidney cells by indirect immunofluorescence
Specificity	Purified virions
Shelf Life	Two years from date of manufacture