

Product Datasheet

Anti-Parainfluenza virus type 3 V1201

Product Name	Anti-Parainfluenza virus type 3
Catalog Number	V1201
Species/Host	Goat
Tested Applications	Indirect immunofluorescence
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with type 2. Minimal reactivity with type 1. Some cross-reactivity with bovine para-3 ($\leq 1:400$) & canine para [1,000]. Non-reactive with HEp-2 cells by indirect immunofluorescence
Specificity	All structural antigens
Shelf Life	Six years from date of manufacture