

## Product Datasheet

### Anti-Influenza A virus (matrix) V1321

<b>Product Name</b>	Anti-Influenza A virus (matrix)
<b>Catalog Number</b>	V1321
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	Indirect immunofluorescence
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified IgG fraction prepared by methods developed by Medix Biochemica
<b>Concentration</b>	5 mg/ml
<b>Storage</b>	Short term (<6 months): 2-8°C Long term: -20°C
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥95%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Non-reactive with the M2 matrix protein. Non-reactive with Hep-2 cells by indirect immunofluorescence
<b>Specificity</b>	Purified matrix protein (M1)
<b>Shelf Life</b>	Six years from date of manufacture