

## Product Datasheet

### Anti-Human IL-17A Monoclonal Antibody, Unconjugated Clone B-A50 853.943.020

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Anti-Human IL-17A Monoclonal Antibody, Unconjugated Clone B-A50     |
| <b>Catalog Number</b>      | 853.943.020   |
| <b>Description</b>         | Mouse monoclonal  |
| <b>Species/Host</b>        | Mouse   |
| <b>Conjugation</b>         | Purified  |
| <b>Tested Applications</b> | FC  |
| <b>Immunogen</b>           | Recombinant human IL-17A  |
| <b>Brand</b>               | Diaclone  |
| <b>Form/Appearance</b>     | Phosphate-buffered saline with 1% BSA and 0.09% sodium azide        |
| <b>Storage</b>             | Stable at +2-8°C for 12 months. For longer storage freeze aliquots. |
| <b>Note</b>                | Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood   |
| <b>Isotype</b>             | Mouse IgG2b Kappa light chain                                       |
| <b>Clonality</b>           | Monoclonal  |
| <b>Specificity</b>         | Recognizes both natural and recombinant IL-17A                      |
| <b>Hybridoma</b>           | Myeloma X63/AG.8653 x Balb/c node cells                             |