

## Product Datasheet

### Anti-Chlamydia trachomatis (MOMP) HRP V1624

<b>Product Name</b>	Anti-Chlamydia trachomatis (MOMP) HRP
<b>Catalog Number</b>	V1624
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Short term (<6 months): 2-8°C Long term: -20°C
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥95%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
<b>Cross Reactivity</b>	Non-reactive with C psittacii and C. pneumoniae in MIF. Negative vs HEp-2 cells and egg yolk sac
<b>Specificity</b>	Major Outer Membrane Protein (MOMP)
<b>Shelf Life</b>	Two years from date of manufacture