

Product Datasheet

Anti-Hepatitis B virus V1806

Product Name	Anti-Hepatitis B virus
Catalog Number	V1806
Species/Host	Horse
Tested Applications	ELISA, gel techniques
Brand	ViroStat
Form/Appearance	Immunoaffinity purified specific antibody
Concentration	1 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with human serum proteins
Specificity	HbsAg
Shelf Life	Six years from date of manufacture