

Product Datasheet

Anti-HIV-1 (gp120) V1961

Product Name	Anti-HIV-1 (gp120)
Catalog Number	V1961
Species/Host	Goat
Tested Applications	WB immunoassay
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with human T or B cells. Non-reactive with human or bovine sera
Specificity	gp120 (predominate gp120 band in WB)
Shelf Life	Six years from date of manufacture