

## Product Datasheet

### Anti-Bacillus HRP V3404

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Anti-Bacillus HRP  |
| <b>Catalog Number</b>      | V3404  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Tested Applications</b> | ELISA  |
| <b>Brand</b>               | ViroStat   |
| <b>Form/Appearance</b>     | Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica |
| <b>Concentration</b>       | 1 mg/ml  |
| <b>Storage</b>             | Short term (<6 months): 2-8°C Long term: -20°C   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Purity</b>              | ≥95%   |
| <b>Buffer</b>              | 0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal                       |
| <b>Cross Reactivity</b>    | Antiserum is unabsorbed and is expected to cross react with various Bacillus species   |
| <b>Specificity</b>         | Various Bacillus species   |
| <b>Shelf Life</b>          | Two years from date of manufacture   |