

## Product Datasheet

### Anti-Campylobacter V4501

**Product Name** Anti-Campylobacter

**Catalog Number** V4501

**Species/Host** Rabbit

**Tested Applications** ELISA

**Brand** ViroStat

**Form/Appearance** Purified from antiserum by protein A chromatography

**Concentration** 5 mg/ml

**Storage** Short term (<6 months): 2-8°C Long term: -20°C

**Isotype** IgG

**Clonality** Polyclonal

**Purity** ≥95%

**Buffer** 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

**Cross Reactivity** Antiserum has not been absorbed & may react with related microorganisms

**Specificity** Whole cells

**Shelf Life** Six years from date of manufacture