

Product Datasheet

Anti-Human IL-17F Monoclonal Antibody, Azide Free Clone B-G46 853.990

Product Name Anti-Human IL-17F Monoclonal Antibody, Azide Free Clone B-G46

Catalog Number 853.990

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide Free

Tested Applications ELISA, ELISpot

Immunogen Recombinant human IL-17F

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2b Kappa light chain

Clonality Monoclonal

Specificity Recognizes both recombinant and natural human IL-17F

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells