

Product Datasheet

Anti-Giardia lamblia FITC V6223

Product Name	Anti-Giardia lamblia FITC
Catalog Number	V6223
Species/Host	Goat
Tested Applications	ELISA, IFA, Indirect immunofluorescence
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with a rod-shaped microorganism found in rodent feces
Specificity	Intact cysts
Shelf Life	Two years from date of manufacture