

Product Datasheet

Anti-Candida albicans V6411

Product Name	Anti-Candida albicans
Catalog Number	V6411
Species/Host	Rabbit
Tested Applications	ELISA, IHC, IFA, Double-diffusion, CIE
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with human serum, urine, & spinal fluid. Antiserum has not been absorbed and does cross react with other yeasts
Specificity	Recognizes numerous proteins in a soluble C. albicans extract (IEP)
Shelf Life	Six years from date of manufacture