

## Product Datasheet

### Anti-Human CD1b/c Monoclonal Antibody, Azide Free Clone B-B5 854.000

**Product Name** Anti-Human CD1b/c Monoclonal Antibody, Azide Free Clone B-B5

**Catalog Number** 854.000

**Description** Mouse monoclonal

**Species/Host** Mouse

**Conjugation** Azide free

**Tested Applications** FC

**Immunogen** Thymus cells and Jurkat cell line

**Brand** Diaclone

**Form/Appearance** Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

**Storage** Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

**Isotype** Mouse IgG1 Kappa light chain

**Clonality** Monoclonal

**Specificity** Recognises a 43-45 kDa protein

**Hybridoma** Myeloma X63/AG.8653 x Balb/c spleen cells