

Product Datasheet

Anti-Mycoplasma pneumoniae HRP V6514

| | |
|----------------------------|--|
| Product Name | Anti-Mycoplasma pneumoniae HRP |
| Catalog Number | V6514 |
| Species/Host | Rabbit |
| Tested Applications | ELISA |
| Brand | ViroStat |
| Form/Appearance | Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica |
| Concentration | 1 mg/ml |
| Storage | Short term (<6 months): 2-8°C Long term: -20°C |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal |
| Cross Reactivity | Non-reactive with human and bovine sera. Antiserum has not been absorbed and may react with related microorganisms |
| Specificity | All antigens |
| Shelf Life | Two years from date of manufacture |