

Product Datasheet

Anti-Staphylococcus aureus HRP V6884

| | |
|----------------------------|--|
| Product Name | Anti-Staphylococcus aureus HRP |
| Catalog Number | V6884 |
| Species/Host | Rabbit |
| Tested Applications | ELISA |
| Brand | ViroStat |
| Form/Appearance | Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica |
| Concentration | 1 mg/ml |
| Storage | Short term (<6 months): 2-8°C Long term: -20°C |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal |
| Cross Reactivity | Antiserum is not absorbed for related microorganisms |
| Specificity | Soluble and structural antigens of the whole bacterium |
| Shelf Life | Two years from date of manufacture |