

Product Datasheet

Anti-Staphylococcus aureus Biotin V6887

Product Name	Anti-Staphylococcus aureus Biotin
Catalog Number	V6887
Species/Host	Rabbit
Tested Applications	ELISA, Indirect immunofluorescence
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is not absorbed for related microorganisms
Specificity	Soluble and structural antigens of the whole bacterium
Shelf Life	Two years from date of manufacture