

Product Datasheet

Anti-Klebsiella FITC V6893

Product Name	Anti-Klebsiella FITC
Catalog Number	V6893
Species/Host	Rabbit
Tested Applications	IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is unabsorbed and may react with related Enterobacteriaceae
Specificity	Polyvalent for multiple Klebsiella antigens
Shelf Life	Two years from date of manufacture