

Product Datasheet

Anti-Herpes simplex virus, type 2 gG) V0129

Product Name	Anti-Herpes simplex virus, type 2 gG)
Catalog Number	V0129
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HSV 2. Non-reactive with HSV1, CMV, EBV, & VZV
Specificity	Glycoprotein G
Shelf Life	Four years from date of manufacture