

Product Datasheet

Anti-Streptococcus pneumoniae V0471

Product Name	Anti-Streptococcus pneumoniae
Catalog Number	V0471
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG3
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Streptococcus pneumoniae common C-polysaccharide (phosphorylcholine epitope)
Shelf Life	Four years from date of manufacture