

## Product Datasheet

### Anti-Rotavirus V0521

<b>Product Name</b>	Anti-Rotavirus
<b>Catalog Number</b>	V0521
<b>Tested Applications</b>	LF, ELISA, WB, Immunostaining cells
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Non-reactive with Influenza A, Influenza B, RSV, Para 1-2-3, Adenovirus, M. pneumoniae, H. pylori, & mammalian cells. Additional cross reactivity data is available for this clone upon request.
<b>Specificity</b>	Recognizes the P41 capsid protein of monkey, porcine, & human isolates
<b>Shelf Life</b>	Four years from date of manufacture