

## Product Datasheet

### Anti-Rotavirus V0531

<b>Product Name</b>	Anti-Rotavirus
<b>Catalog Number</b>	V0531
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with intact virus strains RRV, WA, SA11, & bovine. Non-reactive with Influenza A, Influenza B, RSV, Para 1-2-3, Adenovirus, M. pneumoniae, H. pylori, & mammalian cells
<b>Specificity</b>	Intact virus strains RRV, WA, SA11, & bovine
<b>Matched Pairs</b>	ELISA pairs (capture-conjugate): V0531-V0541, V0531-V0581
<b>Shelf Life</b>	Four years from date of manufacture