

Product Datasheet

Anti-Respiratory Syncytial virus (fusion) V0631

Product Name	Anti-Respiratory Syncytial virus (fusion)
Catalog Number	V0631
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Influenza A & B, Adenovirus, Parainfluenza virus 1-3, Measles, & Mumps
Specificity	Fusion protein of RSV, types A & B
Matched Pairs	ELISA pairs (capture-conjugate): V0631-V0635, V0631-V0651
Shelf Life	Four years from date of manufacture