

## Product Datasheet

### Anti-Respiratory Syncytial virus (fusion) V0635

<b>Product Name</b>	Anti-Respiratory Syncytial virus (fusion)
<b>Catalog Number</b>	V0635
<b>Tested Applications</b>	ELISA, IFA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Non-reactive with Influenza A & B, Adenovirus, Parainfluenza virus 1-3, Measles, & Mumps
<b>Specificity</b>	Fusion protein of RSV, types A & B
<b>Matched Pairs</b>	ELISA pairs (capture-conjugate): V0631-V0635, V0651-V0635, V0641-V0635, V0635-V0651
<b>Shelf Life</b>	Four years from date of manufacture