

Product Datasheet

Anti-Shigella flexneri (2a) V0944

| | |
|----------------------------|---|
| Product Name | Anti-Shigella flexneri (2a) |
| Catalog Number | V0944 |
| Tested Applications | ELISA (whole cells & LPS), IFA (whole cells) |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥ 90% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Non-reactive with Shigella flexneri 1a, 3a, 4a, 5a, 6, Salmonellae (mix), & E. coli (mix) |
| Specificity | LPS of Shigella flexneri 2a |
| Shelf Life | Four years from date of manufacture |