

Product Datasheet

Anti-Parainfluenza virus type 2 V1241

Product Name Anti-Parainfluenza virus type 2

Catalog Number V1241

Tested Applications IFA

Brand ViroStat

Form/Appearance Purified from mouse ascites by protein A chromatography

Concentration ≥ 1 mg/ml

Storage 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.

Isotype IgG2a

Clonality Monoclonal

Purity $\geq 90\%$

Buffer 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

Cross Reactivity Non-reactive with types 1 & 3 by indirect IFA

Specificity Parainfluenza virus type 2

Shelf Life Four years from date of manufacture