

Product Datasheet

Anti-Parainfluenza virus type 3 V1256

Product Name	Anti-Parainfluenza virus type 3
Catalog Number	V1256
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with Bovine PI3.
Specificity	Hemagglutinin of Parainfluenza virus, type 3
Matched Pairs	Pairs in ELISA with V1258
Shelf Life	Four years from date of manufacture