

Product Datasheet

Anti-Adenovirus V1431

Product Name	Anti-Adenovirus
Catalog Number	V1431
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with 34 serotypes of Adenovirus including types 40 & 41. Non-reactive with Influenza A, Influenza B, RSV, Para 1, 2, & 3, Mycoplasma pneumoniae, H. pylori and mammalian cells. Additional cross reactivity data is available upon request.
Specificity	Adenovirus hexon group antigen
Matched Pairs	ELISA pairs (capture-conjugate) : V1431-V1411, V1431-V1461, V1431-V1481, V1411- V1431, V1421-V1461, V1421-V1411
Shelf Life	Four years from date of manufacture