

Product Datasheet

Anti-Chlamydia trachomatis (LPS) V1644

Product Name	Anti-Chlamydia trachomatis (LPS)
Catalog Number	V1644
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with <i>C. trachomatis</i> , <i>C. pneumoniae</i> , <i>C. psittaci</i>
Specificity	LPS
Matched Pairs	ELISA pairs (capture-conjugate): V1645-V1644, V1644-V1644, V1646-V1644, V1647-V1644
Shelf Life	Four years from date of manufacture