

Product Datasheet

Anti-Chlamydia trachomatis (MOMP) V1655

Product Name	Anti-Chlamydia trachomatis (MOMP)
Catalog Number	V1655
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Chlamydia pneumoniae and Feline Chlamydia
Specificity	Major outer membrane protein (MOMP)
Shelf Life	Four years from date of manufacture