

Product Datasheet

Anti-Chlamydia (LPS) V1681

Product Name	Anti-Chlamydia (LPS)
Catalog Number	V1681
Tested Applications	ELISA, IHC, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with <i>C. trachomatis</i> , <i>C. pneumoniae</i> , & <i>C. psittaci</i>
Specificity	LPS
Shelf Life	Four years from date of manufacture